

Antibodies Product Datasheet

Product Name:	SLC22A15 Antibody
Catalog Number:	C50927
Reactivity:	Human
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit, IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SLC22A15
Concentration:	1 mg/ml
Gene Name:	SLC22A15 FLIPT1 UNQ9429/PRO34686
Protein Name:	Solute carrier family 22 member 15 (Fly-like putative transporter 1) (FlIPT1)
Other Name:	
Human Swiss Prot No.:	Q8IZD6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SLC22A15 at Human
Storage and Stability:	Store at -20° C/1 year