

Antibodies Product Datasheet

Product Name:	NPA1 Antibody
Catalog Number:	C17153
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human URB1. AA range:2013-2062
Concentration:	1 mg/ml
Gene Name:	URB1
Protein Name:	Nucleolar pre-ribosomal-associated protein 1
Other Name:	URB1; C21orf108; KIAA0539; NOP254; NPA1; Nucleolar pre-ribosomal-associated protein 1; Nucleolar protein 254 kDa; URB1 ribosome biogenesis 1 homolog
Human Swiss Prot No.:	O60287
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NPA1 Polyclonal Antibody detects endogenous levels of NPA1 protein.
Storage and Stability:	-20°C/1 year