

## Antibodies Product Datasheet

<b>Product Name:</b>	Parathyroid Hormone (PTH) (PT0210) mouse mAb
<b>Catalog Number:</b>	N1803
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	IHC-p,IF(paraffin section)
<b>Dilutions:</b>	IHC-p 1:100-500
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse/IgG1, Kappa
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human Parathyroid Hormone (PTH)
<b>Concentration:</b>	
<b>Gene Name:</b>	PTH
<b>Protein Name:</b>	Parathyroid hormone (PTH) (Parathormone) (Parathyrin)
<b>Other Name:</b>	
<b>Human Swiss Prot No.:</b>	P01270
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of human Parathyroid Hormone (PTH). Heat-induced epitope retrieval (HI-ER) TRIS-EDTA of pH8.0 was highly recommended as antigen repair method in paraffin section
<b>Storage and Stability:</b>	-20°C/1 year