

Antibodies Product Datasheet

Product Name:	HORMAD1 Antibody
Catalog Number:	C11285
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HORMAD1. AA range:321-370
Concentration:	1 mg/ml
Gene Name:	HORMAD1
Protein Name:	HORMA domain-containing protein 1
Other Name:	HORMAD1; NOHMA; HORMA domain-containing protein 1; Cancer/testis antigen 46; CT46; Newborn ovary HORMA protein Q86X24
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HORMAD1 Polyclonal Antibody detects endogenous levels of HORMAD1 protein.
Storage and Stability:	-20°C/1 year