

Antibodies Product Datasheet

Product Name:	TUSC2 Antibody
Catalog Number:	C10196
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TUSC2. AA range:41-90
Concentration:	1 mg/ml
Gene Name:	TUSC2
Protein Name:	Tumor suppressor candidate 2
Other Name:	TUSC2; C3orf11; FUS1; LGCC; PDAP2; Tumor suppressor candidate 2; Fusion 1 protein; Fus-1 protein; PDGFA-associated protein 2
Human Swiss Prot No.:	O75896
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TUSC2 Polyclonal Antibody detects endogenous levels of TUSC2 protein.
Storage and Stability:	-20°C/1 year