

Antibodies Product Datasheet

Product Name:	GDPD5 Antibody
Catalog Number:	C50126
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GDPD5
Concentration:	1 mg/ml
Gene Name:	GDPD5 GDE2 PP6037 PP9363 UNQ1850/PRO3580
Protein Name:	Glycerophosphodiester phosphodiesterase domain-containing protein 5 (EC 3.1.-.-) (Glycerophosphodiester phosphodiesterase 2)
Other Name:	
Human Swiss Prot No.:	Q8WTR4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of GDPD5 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year