

## Antibodies Product Datasheet

<b>Product Name:</b>	OSR1 (Phospho-Thr185) Antibody
<b>Catalog Number:</b>	A8154
<b>Reactivity:</b>	Human, Mouse
<b>Applications:</b>	WB, IHC-p, IF(paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	The antiserum was produced against synthesized peptide derived from human OSR1 around the phosphorylation site of Thr185. AA range:151-200
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	OXSR1
<b>Protein Name:</b>	Serine/threonine-protein kinase OSR1
<b>Other Name:</b>	OXSR1; KIAA1101; OSR1; Serine/threonine-protein kinase OSR1; Oxidative stress-responsive 1 protein
<b>Human Swiss Prot No.:</b>	O95747
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	Phospho-OSR1 (T185) Polyclonal Antibody detects endogenous levels of OSR1 protein only when phosphorylated at T185.
<b>Storage and Stability:</b>	-20°C/1 year