
Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | CD2 Monoclonal Antibody(3F3) |
| Catalog Number: | N1334 |
| Reactivity: | Human,Mouse,Rat |
| Applications: | IHC-p,IF(paraffin section) |
| Dilutions: | IHC: 1:200 |
| Modification: | Unmodified/ Total |
| Source: | Monoclonal Mouse |
| Purification: | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |
| Immunogen: | Synthetic Peptide of CD2 |
| Concentration: | |
| Gene Name: | CD2 |
| Protein Name: | T-cell surface antigen CD2 |
| Other Name: | CD2; SRBC; T-cell surface antigen CD2; Erythrocyte receptor; LFA-2; LFA-3 receptor; Rosette receptor; T-cell surface antigen T11/Leu-5; CD2 |
| Human Swiss Prot No.: | P06729 |
| Formulation: | PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol. |
| Specificity: | The antibody detects endogenous CD2 proteins. |
| Storage and Stability: | -20°C/1 year |