
Antibodies Product Datasheet

Product Name:	c-Fos Antibody
Catalog Number:	B8212
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FOS. AA range:15-64
Concentration:	1 mg/ml
Gene Name:	FOS
Protein Name:	Proto-oncogene c-Fos
Other Name:	FOS; G0S7; Proto-oncogene c-Fos; Cellular oncogene fos; G0/G1 switch regulatory protein 7
Human Swiss Prot No.:	P01100
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	c-Fos Polyclonal Antibody detects endogenous levels of c-Fos protein.
Storage and Stability:	-20°C/1 year