

Antibodies Product Datasheet

Product Name:	STK3 Antibody
Catalog Number:	C46486
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000 ELISA 1:5000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human protein . at AA range: 330-410
Concentration:	1 mg/ml
Gene Name:	STK3 KRS1 MST2
Protein Name:	Serine/threonine-protein kinase 3 (EC 2.7.11.1) (Mammalian STE20-like protein kinase 2) (MST-2) (STE20-like kinase MST2) (Serine/threonine-protein kinase Krs-1) [Cleaved into: Serine/threonine-protein kinase 3 36kDa subunit (MST2/N); Serine/threonine-prot
Other Name:	
Human Swiss Prot No.:	Q13188
Formulation:	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Specificity:	STK3 Polyclonal Antibody detects endogenous levels of protein.
Storage and Stability:	-20°C/1 year