
Antibodies Product Datasheet

| | |
|-------------------------------|--|
| Product Name: | Rabbit Anti-Goat IgG (H&L) - Cy5.5 |
| Catalog Number: | SA1231 |
| Reactivity: | Goat |
| Conjugate: | Cy5.5 |
| Applications: | ELISA, IHC, IF, FCM |
| Dilutions: | IHC (1/200 - 1/1000), IF (1/200 - 1/1000), FCM (1/2000 - 1/4000), ELISA (Use at an assay dependent concentration) |
| Modification: | Not Modification Specific |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Concentration: | |
| Other Name: | |
| Formulation: | Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% BSA, 50% glycerol, 0.02% Sodium Azide |
| Storage and Stability: | -20°C/1 year |