

Antibodies Product Datasheet

Product Name:	Laminin γ -2 Antibody
Catalog Number:	C30224
Reactivity:	Human, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human LAMC2. AA range:1021-1070
Concentration:	1 mg/ml
Gene Name:	LAMC2
Protein Name:	Laminin subunit gamma-2
Other Name:	LAMC2; LAMB2T; LAMNB2; Laminin subunit gamma-2; Cell-scattering factor 140 kDa subunit; CSF 140 kDa subunit; Epiligrin subunit gamma; Kalinin subunit gamma; Kalinin/nicein/epiligrin 100 kDa subunitLadsin 140 kDa subunit; Laminin B2t chain; Laminin-5 subun
Human Swiss Prot No.:	Q13753
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Laminin γ -2 Polyclonal Antibody detects endogenous levels of Laminin γ -2 protein.
Storage and Stability:	-20°C/1 year