

Antibodies Product Datasheet

Product Name:	MARK3 Antibody
Catalog Number:	R12-2998
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC, ELISA
Dilutions:	MARK3 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MARK3.
Concentration:	1 mg/ml
Gene Name:	MARK3
Other Name:	C-TAK1; Cdc25C- associated protein kinase 1; MAP/microtubule affinity-regulating kinase 3; Protein kinase STK10; Ser/Thr protein kinase PAR-1; Serine/threonine protein kinase p78; cTAK1
Human Swiss Prot No.:	P27448
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	MARK3 Antibody antibody detects endogenous levels of MARK3.
Storage and Stability:	Stable at -20°C for at least 1 year.