

Antibodies Product Datasheet

Product Name:	TBX10 Antibody
Catalog Number:	C10164
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TBX10. AA range:261-310
Concentration:	1 mg/ml
Gene Name:	TBX10
Protein Name:	T-box transcription factor TBX10
Other Name:	TBX10; TBX7; T-box transcription factor TBX10; T-box protein 10
Human Swiss Prot No.:	O75333
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TBX10 Polyclonal Antibody detects endogenous levels of TBX10 protein.
Storage and Stability:	-20°C/1 year