

Antibodies Product Datasheet

Product Name:	CPSF7 Antibody
Catalog Number:	C49959
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human CPSF7
Concentration:	1 mg/ml
Gene Name:	CPSF7
Protein Name:	Cleavage and polyadenylation specificity factor subunit 7 (Cleavage and polyadenylation specificity factor 59 kDa subunit) (CFIm59) (CPSF 59 kDa subunit) (Pre-mRNA cleavage factor Im 59 kDa subunit)
Other Name:	
Human Swiss Prot No.:	Q8N684
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of CPSF7 at Human, Mouse,Rat
Storage and Stability:	Store at -20° C/1 year