

## Antibodies Product Datasheet

<b>Product Name:</b>	p53 (Acetyl-K305) Antibody
<b>Catalog Number:</b>	D12169
<b>Reactivity:</b>	Human:K305
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	wb dilution 1:1000
<b>Modification:</b>	Acetyl Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Synthesized Acetyl peptide derived from human p53. at AA range: K305
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	TP53 P53
<b>Protein Name:</b>	p53
<b>Other Name:</b>	Cellular tumor antigen p53 (Antigen NY-CO-13) (Phosphoprotein p53) (Tumor suppressor p53)
<b>Human Swiss Prot No.:</b>	P04637
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of p53 (Acetyl-K305), It doesn't react with total protein.
<b>Storage and Stability:</b>	-20°C/1 year