

Antibodies Product Datasheet

Product Name:	Siah-1/2 Antibody
Catalog Number:	C11300
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SIAH1. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	SIAH1/SIAH2
Protein Name:	E3 ubiquitin-protein ligase SIAH1/2
Other Name:	SIAH1; HUMSIAH; E3 ubiquitin-protein ligase SIAH1; Seven in absentia homolog 1; Siah-1; Siah-1a; SIAH2; E3 ubiquitin-protein ligase SIAH2; Seven in absentia homolog 2; Siah-2; hSiah2
Human Swiss Prot No.:	Q8IUQ4/O43255
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Siah-1/2 Polyclonal Antibody detects endogenous levels of Siah-1/2 protein.
Storage and Stability:	-20°C/1 year