

---

## Antibodies Product Datasheet

<b>Product Name:</b>	SRT-Tag Monoclonal Antibody(4H9), HRP Conjugated
<b>Catalog Number:</b>	N1201
<b>Reactivity:</b>	Species independent
<b>Applications:</b>	WB
<b>Dilutions:</b>	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:5000.
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse IgG
<b>Purification:</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	
<b>Concentration:</b>	1mg/ml
<b>Gene Name:</b>	
<b>Protein Name:</b>	
<b>Other Name:</b>	SRT-Tag
<b>Human Swiss Prot No.:</b>	
<b>Formulation:</b>	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
<b>Specificity:</b>	SRT-Tag Monoclonal Antibody(4H9) HRP Conjugated, specially designed for your Western blot analysis.
<b>Storage and Stability:</b>	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.