

Antibodies Product Datasheet

Product Name:	SLC10A4 Antibody
Catalog Number:	C18824
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from SLC10A4 . at AA range: 340-420
Concentration:	1 mg/ml
Gene Name:	SLC10A4
Protein Name:	Sodium/bile acid cotransporter 4
Other Name:	SLC10A4; Sodium/bile acid cotransporter 4; Na(+)/bile acid cotransporter 4; Solute carrier family 10 member 4
Human Swiss Prot No.:	Q96EP9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SLC10A4 Polyclonal Antibody detects endogenous levels of SLC10A4 protein.
Storage and Stability:	-20°C/1 year