

Antibodies Product Datasheet

Product Name:	FoxA2 Monoclonal Antibody
Catalog Number:	N0280
Reactivity:	Human
Applications:	WB,IHC-p,IF(paraffin section),ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of FoxA2 expressed in E. Coli.
Concentration:	
Gene Name:	FOXA2
Protein Name:	Hepatocyte nuclear factor 3-beta
Other Name:	FOXA2; HNF3B; TCF3B; Hepatocyte nuclear factor 3-beta; HNF-3-beta; HNF-3B; Forkhead box protein A2; Transcription factor 3B; TCF-3B
Human Swiss Prot No.:	Q9Y261
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	FoxA2 Monoclonal Antibody detects endogenous levels of FoxA2 protein.
Storage and Stability:	-20°C/1 year