

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | NEDD4L Antibody |
| Catalog Number: | C44963 |
| Reactivity: | Human, Mouse |
| Applications: | WB |
| Dilutions: | WB 1:500-2000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Immunogen: | Synthesized peptide derived from human NEDD4L |
| Concentration: | 1 mg/ml |
| Gene Name: | NEDD4L KIAA0439 NEDL3 |
| Protein Name: | NEDD4L |
| Other Name: | E3 ubiquitin-protein ligase NEDD4-like (EC 6.3.2.-) (NEDD4.2) (Nedd4-2) |
| Human Swiss Prot No.: | Q96PU5 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.08% sodium azide. |
| Specificity: | This antibody detects endogenous levels of NEDD4L at Human/Mouse |
| Storage and Stability: | -20°C/1 year |