

Antibodies Product Datasheet

Product Name:	p53 (Tri Methyl-Lys370) Antibody
Catalog Number:	D12867
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Meythl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human p53 (Tri Methyl Lys370)
Concentration:	1 mg/ml
Gene Name:	TP53 P53
Protein Name:	p53 (Tri Methyl Lys370)
Other Name:	Cellular tumor antigen p53 (Antigen NY-CO-13) (Phosphoprotein p53) (Tumor suppressor p53)
Human Swiss Prot No.:	P04637
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat p53 (Tri Methyl Lys370)
Storage and Stability:	-20°C/1 year