
Antibodies Product Datasheet

Product Name:	PD1 Nano antibody
Catalog Number:	NA0001
Reactivity:	Human, Mouse
Applications:	FCM, ELISA
Dilutions:	ELISA 1:10000-20000
Modification:	Not Modification Specific
Source:	Nano-antibody with his-tag, Camel
Purification:	The antibody was affinity-purified nickel column using his-tag in antibody
Immunogen:	Recombinant protein of full length human PD1
Concentration:	
Gene Name:	PDCD1 PD1
Protein Name:	Programmed cell death protein 1 (Protein PD-1) (hPD-1) (CD antigen CD279)
Other Name:	
Human Swiss Prot No.:	Q15116
Formulation:	Liquid in PBS
Specificity:	This antibody detects endogenous levels of PD1
Storage and Stability:	-20°C/1 year