

Antibodies Product Datasheet

Product Name:	CD19 (Phospho-Tyr531) Antibody
Catalog Number:	A0412
Reactivity:	Human, Mouse, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CD19 around the phosphorylation site of Tyr531. AA range:501-550
Concentration:	1 mg/ml
Gene Name:	CD19
Protein Name:	B-lymphocyte antigen CD19
Other Name:	CD19; B-lymphocyte antigen CD19; B-lymphocyte surface antigen B4; Differentiation antigen CD19; T-cell surface antigen Leu-12; CD antigen CD19
Human Swiss Prot No.:	P15391
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-CD19 (Y531) Polyclonal Antibody detects endogenous levels of CD19 protein only when phosphorylated at Y531.
Storage and Stability:	-20°C/1 year