

---

## Antibodies Product Datasheet

<b>Product Name:</b>	TTR mouse Monoclonal Antibody(5A11)
<b>Catalog Number:</b>	N1677
<b>Reactivity:</b>	Human
<b>Applications:</b>	WB,IHC-p,IF(paraffin section)
<b>Dilutions:</b>	WB 1:500-2000,IHC-p 1:50-300
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Immunogen:</b>	Recombinant Protein of TTR
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	TTR PALB
<b>Protein Name:</b>	Transthyretin (ATTR) (Prealbumin) (TBPA)
<b>Other Name:</b>	Transthyretin (ATTR) (Prealbumin) (TBPA)
<b>Human Swiss Prot No.:</b>	P02766
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	The antibody detects endogenous TTR protein
<b>Storage and Stability:</b>	-20°C/1 year