

Antibodies Product Datasheet

| | |
|-------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name: | NOTCH3 (Cleaved-Asp1628) Antibody |
| Catalog Number: | L0443 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB, ELISA |
| Dilutions: | WB 1:1000-2000 ELISA 1:5000-20000 |
| Modification: | Cleaved Specific |
| Source: | Polyclonal Rabbit |
| Purification: | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Immunogen: | Synthesized peptide derived from human NOTCH3 (Cleaved-Asp1628) |
| Concentration: | 1 mg/ml |
| Gene Name: | NOTCH3 |
| Protein Name: | NOTCH3 (Cleaved-Asp1628) |
| Other Name: | Neurogenic locus notch homolog protein 3 (Notch 3) [Cleaved into: Notch 3 extracellular truncation; Notch 3 intracellular domain] |
| Human Swiss Prot No.: | Q9UM47 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | This antibody detects endogenous levels of Human, Mouse, Rat NOTCH3 (Cleaved-Asp1628, protein was cleaved amino acid sequence between 1628-1629) |
| Storage and Stability: | -20°C/1 year |