

Antibodies Product Datasheet

Product Name:	EPB41 (Phospho-Tyr660/418) Antibody
Catalog Number:	P12-1101
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	EPB41 (Phospho-Tyr660/418) Antibody
Modification:	Phospho Specific
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using phospho peptide. The antibody against non-phospho peptide was removed by chromatography using corresponding non-phospho peptide.
Immunogen:	The antiserum was produced against synthesized phospho-peptide derived from human EPB41 around the phosphorylation site of Tyr660/418.
Concentration:	1 mg/ml
Gene Name:	EPB41
Other Name:	4.1R; 41; band 4.1; E41P; EPB4.1; erythrocyte membrane protein band 4.1; P4.1; protein 4.1
Human Swiss Prot No.:	P11171
Formulation:	Rabbit IgG in phosphate buffered saline (Phospho-without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	EPB41 (Phospho-Tyr660/Tyr418) Antibody detects endogenous levels of EPB41 only when phosphorylated at Tyr660/Tyr418.
Storage and Stability:	Stable at -20°C for at least 1 year.