
Antibodies Product Datasheet

Product Name:	AKAP14 Antibody
Catalog Number:	R12-2470
Reactivity:	Human
Applications:	WB, ELISA
Dilutions:	AKAP14 Antibody
Modification:	Unmodified/Total
Source:	Rabbit Polyclonal
Purification:	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AKAP14.
Concentration:	1 mg/ml
Gene Name:	AKAP14
Other Name:	A-kinase anchor protein 14; A-kinase anchor protein 14; A-kinase anchor protein 28 kDa; AKAP14; AKAP28
Human Swiss Prot No.:	Q86UN6
Formulation:	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Specificity:	AKAP14 Antibody antibody detects endogenous levels of AKAP14.
Storage and Stability:	Stable at -20°C for at least 1 year.