

Antibodies Product Datasheet

Product Name:	ENTPD4 Antibody
Catalog Number:	C50915
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human ENTPD4
Concentration:	1 mg/ml
Gene Name:	ENTPD4 KIAA0392 LALP70 LYSAL1
Protein Name:	Ectonucleoside triphosphate diphosphohydrolase 4 (NTPDase 4) (EC 3.6.1.6) (Lysosomal apyrase-like protein of 70 kDa) (Uridine-diphosphatase) (UDPase)
Other Name:	
Human Swiss Prot No.:	Q9Y227
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of ENTPD4 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year