

Antibodies Product Datasheet

Product Name:	Nkx-6.3 Antibody
Catalog Number:	C10009
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NKX6.3. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	NKX6-3
Protein Name:	Homeobox protein Nkx-6.3
Other Name:	NKX6-3; Homeobox protein Nkx-6.3
Human Swiss Prot No.:	A6NJ46
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Nkx-6.3 Polyclonal Antibody detects endogenous levels of Nkx-6.3 protein.
Storage and Stability:	-20°C/1 year