

Antibodies Product Datasheet

Product Name:	TSC-22 Antibody
Catalog Number:	C11099
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TSC22D1. AA range:71-120
Concentration:	1 mg/ml
Gene Name:	TSC22D1
Protein Name:	TSC22 domain family protein 1
Other Name:	TSC22D1; KIAA1994; TGFB114; TSC22; hucep-2; TSC22 domain family protein 1; Cerebral protein 2; Regulatory protein TSC-22; TGFB-stimulated clone 22 homolog; Transforming growth factor beta-1-induced transcript 4 protein
Human Swiss Prot No.:	Q15714
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TSC-22 Polyclonal Antibody detects endogenous levels of TSC-22 protein.
Storage and Stability:	-20°C/1 year