

Antibodies Product Datasheet

Product Name:	VHL Antibody
Catalog Number:	C30872
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p 1:50-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the N-terminal region of human VHL. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	VHL
Protein Name:	Von Hippel-Lindau disease tumor suppressor (Protein G7) (pVHL)
Other Name:	Von Hippel-Lindau disease tumor suppressor (Protein G7) (pVHL)
Human Swiss Prot No.:	P40337
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous VHL
Storage and Stability:	-20°C/1 year