

Antibodies Product Datasheet

Product Name:	XRCC4 (Phospho-Ser260) Antibody
Catalog Number:	A8452
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human XRCC4 (Phospho-Ser260)
Concentration:	1 mg/ml
Gene Name:	XRCC4
Protein Name:	XRCC4 (Phospho-Ser260)
Other Name:	DNA repair protein XRCC4 (X-ray repair cross-complementing protein 4)
Human Swiss Prot No.:	Q13426
Formulation:	Liquid in PBS containing 50% glycerol, and 0.151% sodium azide.
Specificity:	This antibody detects endogenous levels of XRCC4 (Phospho-Ser260) at Human, Mouse,Rat
Storage and Stability:	-20°C/1 year