

## Antibodies Product Datasheet

<b>Product Name:</b>	MMP9 (Cleaved-Phe107) Antibody
<b>Catalog Number:</b>	L0470
<b>Reactivity:</b>	Human, Rat, Mouse
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB 1:1000-2000 ELISA 1:5000-20000
<b>Modification:</b>	Cleaved Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized peptide derived from human MMP9 (Cleaved-Phe107)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	MMP9 CLG4B
<b>Protein Name:</b>	MMP9 (Cleaved-Phe107)
<b>Other Name:</b>	Matrix metalloproteinase-9 (MMP-9) (EC 3.4.24.35) (92 kDa gelatinase) (92 kDa type IV collagenase) (Gelatinase B) (GELB) [Cleaved into: 67 kDa matrix metalloproteinase-9; 82 kDa matrix metalloproteinase-9]
<b>Human Swiss Prot No.:</b>	P14780
<b>Formulation:</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Specificity:</b>	This antibody detects endogenous levels of Human MMP9 (Cleaved-Phe107, protein was cleaved amino acid sequence between 106-107 )
<b>Storage and Stability:</b>	-20°C/1 year