

---

## Antibodies Product Datasheet

<b>Product Name:</b>	BMPR-II Monoclonal Antibody
<b>Catalog Number:</b>	N0073
<b>Reactivity:</b>	Human, Mouse, Rat, Monkey
<b>Applications:</b>	WB, IHC-p, IF (paraffin section), ELISA
<b>Dilutions:</b>	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
<b>Modification:</b>	Unmodified/ Total
<b>Source:</b>	Monoclonal Mouse
<b>Purification:</b>	Affinity purification
<b>Immunogen:</b>	Purified recombinant fragment of human BMPR-II expressed in E. Coli.
<b>Concentration:</b>	
<b>Gene Name:</b>	BMPR2
<b>Protein Name:</b>	Bone morphogenetic protein receptor type-2
<b>Other Name:</b>	BMPR2; PPH1; Bone morphogenetic protein receptor type-2; BMP type-2 receptor; BMPR-2; Bone morphogenetic protein receptor type II; BMP type II receptor; BMPR-II
<b>Human Swiss Prot No.:</b>	Q13873
<b>Formulation:</b>	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
<b>Specificity:</b>	BMPR-II Monoclonal Antibody detects endogenous levels of BMPR-II protein.
<b>Storage and Stability:</b>	-20°C/1 year