

Antibodies Product Datasheet

Product Name:	SPACA3 Antibody
Catalog Number:	C50083
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	WB 1:500-2000
Modification:	Unmodified/Total
Source:	Polyclonal, Rabbit,IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human SPACA3
Concentration:	1 mg/ml
Gene Name:	SPACA3 LYC3 LYZL3 SLLP1 SPRASA UNQ424/PRO862
Protein Name:	Sperm acrosome membrane-associated protein 3 (Cancer/testis antigen 54) (CT54) (Lysozyme-like acrosomal sperm-specific secretory protein ALLP-17) (Lysozyme-like protein 3) (Sperm lysozyme-like protein 1) (Sperm protein reactive with antisperm antibodies)
Other Name:	
Human Swiss Prot No.:	Q8IXA5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SPACA3 at Human, Mouse
Storage and Stability:	Store at -20° C/1 year