

Antibodies Product Datasheet

Product Name:	SR-2C Antibody
Catalog Number:	C12015
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human 5-HT-2C. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	HTR2C
Protein Name:	5-hydroxytryptamine receptor 2C
Other Name:	HTR2C; HTR1C; 5-hydroxytryptamine receptor 2C; 5-HT-2C; 5-HT2C; 5-HT2C; 5-hydroxytryptamine receptor 1C; 5-HT-1C; 5-HT1C; Serotonin receptor 2C
Human Swiss Prot No.:	P28335
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	SR-2C Polyclonal Antibody detects endogenous levels of SR-2C protein.
Storage and Stability:	-20°C/1 year