

## Antibodies Product Datasheet

<b>Product Name:</b>	MER/TYRO3 (Phospho-Tyr753/Tyr685) Antibody
<b>Catalog Number:</b>	A9045
<b>Reactivity:</b>	Human:Y753/Y685, Mouse:Y748/Y675, Rat:Y748/Y675
<b>Applications:</b>	WB, ELISA
<b>Dilutions:</b>	WB: 1:500-10000 ELISA: 1:10000
<b>Modification:</b>	Phospho Specific
<b>Source:</b>	Polyclonal Rabbit
<b>Purification:</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Immunogen:</b>	Synthesized phospho-peptide around the phosphorylation site of human MER/TYRO3 (Phospho-Tyr753/Tyr685)
<b>Concentration:</b>	1 mg/ml
<b>Gene Name:</b>	MERTK MER
<b>Protein Name:</b>	c-mer proto-oncogene tyrosine kinase
<b>Other Name:</b>	MERTK MER
<b>Human Swiss Prot No.:</b>	Q12866/Q06418
<b>Formulation:</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Specificity:</b>	The antibody detects endogenous MER/TYRO3 (Phospho-Tyr753/Tyr685) protein
<b>Storage and Stability:</b>	-20°C/1 year