
Antibodies Product Datasheet

Product Name:	BPTF Monoclonal Antibody
Catalog Number:	N0079
Reactivity:	Human
Applications:	WB,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Purified recombinant fragment of human BPTF expressed in E. Coli.
Concentration:	
Gene Name:	BPTF
Protein Name:	Nucleosome-remodeling factor subunit BPTF
Other Name:	BPTF; FAC1; FALZ; Nucleosome-remodeling factor subunit BPTF; Bromodomain and PHD finger-containing transcription factor; Fetal Alz-50 clone 1 protein; Fetal Alzheimer antigen
Human Swiss Prot No.:	Q12830
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	BPTF Monoclonal Antibody detects endogenous levels of BPTF protein.
Storage and Stability:	-20°C/1 year