

Antibodies Product Datasheet

Product Name:	ATRIP Antibody
Catalog Number:	B0772
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ATRIP. AA range:34-83
Concentration:	1 mg/ml
Gene Name:	ATRIP
Protein Name:	ATR-interacting protein
Other Name:	ATRIP; AGS1; ATR-interacting protein; ATM and Rad3-related-interacting protein
Human Swiss Prot No.:	Q8WXE1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATRIP Polyclonal Antibody detects endogenous levels of ATRIP protein.
Storage and Stability:	-20°C/1 year