

Antibodies Product Datasheet

Product Name:	HMGN1/2/3/4 (Acetyl-Lys27/K33/K31) Antibody
Catalog Number:	D12312
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:1000-2000 ELISA 1:5000-20000
Modification:	Acetyl Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human HMGN1/2/3/4 (Acetyl Lys27/K33/K31)
Concentration:	1 mg/ml
Gene Name:	HMGN1 HMG14
Protein Name:	HMGN1/2/3/4 (Acetyl Lys27/K33/K31)
Other Name:	Non-histone chromosomal protein HMG-14 (High mobility group nucleosome-binding domain-containing protein 1)
Human Swiss Prot No.:	P05114/P05204/Q15651/O00479
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of Human, Mouse, Rat HMGN1/2/3/4 (Acetyl Lys27/K33/K31)
Storage and Stability:	-20°C/1 year