

Antibodies Product Datasheet

Product Name:	FGFR-4 (Phospho-Tyr642) Antibody
Catalog Number:	A9041
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Phospho Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized phospho-peptide around the phosphorylation site of human FGFR-4 (phospho Tyr642)
Concentration:	1 mg/ml
Gene Name:	FGFR4
Protein Name:	Fibroblast growth factor receptor 4
Other Name:	FGFR4; JTK2; TKF; Fibroblast growth factor receptor 4; FGFR-4; CD antigen CD334
Human Swiss Prot No.:	P22455
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Phospho-FGFR-4 (Y642) Polyclonal Antibody detects endogenous levels of FGFR-4 protein only when phosphorylated at Y642.
Storage and Stability:	-20°C/1 year