

Antibodies Product Datasheet

Product Name:	TM230 Antibody
Catalog Number:	C41697
Reactivity:	Human, Mouse, Rat
Applications:	WB
Dilutions:	WB 1: 500-2000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human TM230
Concentration:	1 mg/ml
Gene Name:	TMEM230 C20orf30 HSPC274 UNQ2432/PRO4992
Protein Name:	TM230
Other Name:	
Human Swiss Prot No.:	Q96A57
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of TM230 at Human/Mouse/Rat
Storage and Stability:	-20°C/1 year