

Antibodies Product Datasheet

Product Name:	Calnexin Monoclonal Antibody
Catalog Number:	N0088
Reactivity:	Human
Applications:	WB,IF/ICC,ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	Affinity purification
Immunogen:	Synthetic peptide corresponding to aa (CEAAEERPWLWVVYILTVAL) of human Calnexin, conjugated to KLH.
Concentration:	
Gene Name:	CANX
Protein Name:	Calnexin
Other Name:	CANX; Calnexin; IP90; Major histocompatibility complex class I antigen-binding protein p88; p90
Human Swiss Prot No.:	P27824
Formulation:	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Specificity:	Calnexin Monoclonal Antibody detects endogenous levels of Calnexin protein.
Storage and Stability:	-20°C/1 year