

Antibodies Product Datasheet

| | |
|-------------------------------|---|
| Product Name: | EIF2S3 Antibody |
| Catalog Number: | C49076 |
| Reactivity: | Human, Mouse, Rat |
| Applications: | WB |
| Dilutions: | WB 1:500-2000 |
| Modification: | Unmodified/Total |
| Source: | Polyclonal, Rabbit,IgG |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Immunogen: | Synthesized peptide derived from human EIF2S3 |
| Concentration: | 1 mg/ml |
| Gene Name: | EIF2S3 EIF2G |
| Protein Name: | Eukaryotic translation initiation factor 2 subunit 3 (Eukaryotic translation initiation factor 2 subunit gamma X) (eIF-2-gamma X) (eIF-2gX) |
| Other Name: | |
| Human Swiss Prot No.: | P41091 |
| Formulation: | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Specificity: | This antibody detects endogenous levels of EIF2S3 at Human, Mouse,Rat |
| Storage and Stability: | Store at -20° C/1 year |