

Antibodies Product Datasheet

Product Name:	E-Cadherin(PT2238) mouse mAb
Catalog Number:	N1841
Reactivity:	Human, Mouse
Applications:	IHC-p,IF(paraffin section)
Dilutions:	IHC-p 1:100-500
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse/IgG1, Kappa
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Immunogen:	Synthesized peptide derived from human E-Cadherin
Concentration:	
Gene Name:	CDH1 CDHE UVO
Protein Name:	Cadherin-1 (CAM 120/80) (Epithelial cadherin) (E-cadherin) (Uvomorulin) (CD antigen CD324) [Cleaved into: E-Cad/CTF1; E-Cad/CTF2; E-Cad/CTF3]
Other Name:	
Human Swiss Prot No.:	P12830
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of human Cadherin-1 (CAM 120/80) (Epithelial cadherin) (E-cadherin) (Uvomorulin) (CD antigen CD324) [Cleaved into: E-Cad/CTF1; E-Cad/CTF2; E-Cad/CTF3]. Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH8.0 was hi
Storage and Stability:	-20°C/1 year