

Antibodies Product Datasheet

Product Name:	Anti-Human Beta-catenin (ABT-Catenin) Mouse Monoclonal Antibody
Catalog Number:	V0154
Clone ID:	ABT-Catenin
Reactivity:	Human, Mouse, Rat
Applications:	IHC, WB
Application Dilutions	IHC 1:100-200, WB 1:500-1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
UniProt:	CTNNB1
Gene Name:	P35222
Other Name:	b-catenin, Beta catenin, Beta-catenin, Cadherin associated protein, Catenin (cadherin associated protein), beta 1, 88 kDa, Catenin beta 1, Catenin beta-1, CATNB, CHBCAT, CTNB1_HUMAN, CTNNB, CTNNB1, DKFZp686D02253, FLJ25606, FLJ37923, OTTHUMP00000162082, OTTHUMP00000165222, OTTHUMP00000165223, OTTHUMP00000209288, OTTHUMP00000209289
Isotype	IgG2b, Kappa
Localization:	Colon/Appendix/Liver
Positive Tissue Control:	Membranous
Immunogen	Synthetic peptide
Purification	Immunogen affinity purified
Formulation:	PBS, pH7.4, 50% glycerol, 0.03%Proclin300
Specificity:	The antibody can specifically recognize human β - Catenin protein. In western blotting of Hela, A431 and HEK293 cell lysates, the antibody labels β -Catenin, no reaction with α or γ -catenin.
Storage and Stability:	-20°C/1 year