

Antibodies Product Datasheet

Product Name:	Caspase-8 Antibody
Catalog Number:	C47571
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section)
Dilutions:	WB: 1:1000-2000 IHC: 1:50-100
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Recombinant Protein of Caspase-8
Concentration:	
Gene Name:	CASP8
Protein Name:	Caspase8
Other Name:	CASP8; MCH5; Caspase-8; CASP-8; Apoptotic cysteine protease; Apoptotic protease Mch-5; CAP4; FADD-homologous ICE/ced-3-like protease; FADD-like ICE; FLICE; ICE-like apoptotic protease 5; MORT1-associated ced-3 homolog; MACH; YM3439
Human Swiss Prot No.:	Q14790
Formulation:	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Specificity:	The antibody detects endogenous Caspase-8 protein.
Storage and Stability:	-20°C/1 year