

Antibodies Product Datasheet

Product Name:	UHRF1(N-term) mouse mAb
Catalog Number:	N0767
Reactivity:	Human, Mouse
Applications:	WB
Dilutions:	wb 1:1000
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human UHRF1 protein fragments expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	uhrf1
Protein Name:	
Other Name:	Ac2-121;AL022808;E3 ubiquitin-protein ligase UHRF1;EC 6.3.2.-;FLJ21925;Hnp95 Huhrf1; HuNp95;ICBP90;Inverted CCAAT box binding protein of 90 kDa;Inverted CCAAT box binding protein, 90-Kd;Inverted CCAAT box-binding protein of 90 kDa;Liver regeneration-relat
Human Swiss Prot No.:	Q96T88
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of UHRF1 and does not cross-react with related proteins
Storage and Stability:	-20°C/1 year