

Antibodies Product Datasheet

Product Name:	SIRT6 mouse mAb
Catalog Number:	N0840
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IF/ICC, IP
Dilutions:	wb 1:500 IF/icc 1:100
Modification:	Unmodified/ Total
Source:	Monoclonal Mouse
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Purified recombinant human SIRT6 protein expressed in E.coli.
Concentration:	1 mg/ml
Gene Name:	sirt6
Protein Name:	
Other Name:	2810449N18Rik; AI043036; Mono ADP ribosyltransferase sirtuin 6; NAD-dependent protein deacetylase sirtuin-6; Regulatory protein SIR2 homolog 6; Regulatory protein SIR2 homolog; SIR2 like 6; SIR2 like protein 6; Sir2 related protein type 6; SIR2-like protein 6; SIR2
Human Swiss Prot No.:	Q8N6T7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of SIRT6 and does not cross-react with related proteins.
Storage and Stability:	-20°C/1 year